



TAPE TEST for pinworms

- Some transparent tape and a slide are needed (the slide can be collected at the CUP).
- Collect the sample in the morning, after having woken up, without cleaning first. Do not apply creams on the anus the previous evening.
- Cut a 5 cm long piece of tape and press the adhesive part on the anus for some seconds.



Press the adhesive part of the tape on the slide, being careful not to leave wrinkles or air bubbles.

It is possible to collect 3 different slides every other day; the slides are to be preserved at room temperature.



**Please note! Do not
use opaque or colored tape!**



Translation by Davide Campari